



## PHYSIOTHERAPY TABLE THOR. REF. 2252A



High performance bed with breathing hole suitable for physiotherapy and osteopathy treatments. It is divided into 8 sections that allow the easy regulation of the table in several positions. The leg supports are independent. It is equipped with 2 motors that control height and backrest. Sturdy metal frame with retractable wheels to facilitate its movement. White PU upholstery of high quality and easy cleaning on high density padding.



## **TECHNICAL DATA**

**DIMENSIONS:** 207\*68\*61-88 cm

**PACKING SIZE:** 186 x 72 x 65 cm

**GROSS WEIGHT / NET WEIGHT:** 86 kg / 69 kg

**VOLTAGE:** 100V-120V/ 220V—240V

**FREQUENCY:** 50/60 HZ

**UPHOLSTERY:** High quality PU

**NUMBER OF MOTORS:** 2

**HANDCONTROLLER:** YES

**BREATHER HOLE:** YES

**ROLL HOLDER:** NO

**ADJUSTABLE ARMREST:** YES (see technical drawing)

**WHEELS:** YES, RETRAC TABLE

**HEATING:** NOT AVAILABLE



**COLORS AVAILABLE:** WHITE, GREEN, BLACK

**MAX. HEIGHT LEGREST EXT.:** 141 cm



**MAXIMUM LOAD:** 175 KG



## BACKREST ADJUSTMENT







## HEADREST ADJUSTMENT





## FOOTREST ADJUSTMENT





## CAUTIONS

- 1) To avoid coming into collision between the different parts, the backrest shouldn't go down at the same time that the headrest when the bed is in the lowest position.
- 2) If the bed is in the lowest position and the backrest is placed horizontally (as in the picture below), the headrest should only go down until  $-22^\circ$ .
- 3) If you are moving the bed into the lowest position, please be careful if the footrest is by  $-10^\circ$  to avoid coming into collision with the floor.

